



PowerTape™

NEW. INNOVATIVE. DIFFERENT.



REGULAR OR BASIC
PREVENTIVE ANKLE USING
2.75" OR 2" POWERFLEX,
1.5" OR 2" POWERTAPE

POWERTAPE
2 - PATENTS PENDING




POST INJURY
TURF TOE/ANKLE COMBINATION USING
1" AND 2" POWERFLEX AND 1" AND 1.5" POWERTAPE



POST-INJURY ELBOW USING
2.75" OR 3" POWERFLEX AND
1.5" OR 2" POWERTAPE



BASIC PREVENTIVE WRIST/HAND/THUMB
USING 1.5" OR 2" POWERFLEX AND
1" AND 1.5" POWERTAPE



POST INJURY TURF TOE
USING 1" AND 2" POWERFLEX
AND 1" OR 1.5" POWERTAPE



SPEED SPATTING USING 2.75" OR 3" POWERFLEX
AND 1.5" OR 2" POWERTAPE

The PowerFlex Taping System
JUST GOT BETTER.



PowerTape™

NEW. INNOVATIVE. DIFFERENT.



Our Philosophy

Taping and wrapping is a skilled art form. Every taping procedure applied by a skilled professional should be a masterpiece...a work of art that an athlete can trust.

Everyday at Andover we are challenged to develop new and innovative products that meet the changing needs of our customers. We strive to deliver products and services that make a difference - we believe bringing the new technology in PowerTape to the sports industry will do just that.

New products and new methods take time to embrace - it requires an open mind and breaking with tradition. Yet, without embracing new products and new technology we can not move forward. We ask that you open your mind to this new taping system and join us as we bring athletic taping to the next level.

Andover is committed to making a difference in the sports medicine industry through new, innovative products and services. We evaluate every idea with an open mind, using imagination and creativity to solve problems and issues that exist in the market today.

PowerTape™

Use in combination with PowerFlex for The Ultimate Taping System.

What is PowerTape?

- PowerTape replaces traditional cloth trainers tape.
- PowerTape is a Cohesive Trainers Tape that sticks only to itself or to PowerFlex.

Why is PowerTape better than traditional cloth tape?

- Creates the STRONGEST combination taping system in existence.
- Contains less than .5% stretch - eliminating the "give" that exists when using traditional trainers tape.
- Remains supportive and stable throughout games and practices.
- Molded training system with better conformability and comfort.
- Holds up against sweat and water.
- Exceptional "lay in".
- Two patents pending.



Andover Coated Products, Inc.

9 Fanaras Drive, Salisbury, MA 01952

Phone: 978.465.0044 • Fax: 978.462.0003 • Toll Free: 800.432.6686

www.andovercoated.com



Proudly Made in U.S.A.

POWERTAPE: 2 - PATENTS PENDING POWERFLEX: U.S. PATENT NO. 5.762.623

Our History

Andover Coated Products has been one of the leading manufacturers of cohesive bandages and tapes for over 25 years and currently supplies the sports medicine, healthcare, and animal health industries with high quality products made in the USA. The company was founded in Massachusetts in 1976, and now distributes its patented products around the world.

In 1995 Andover introduced PowerFlex - a superior cohesive bandage. In 2005 Andover introduced PowerTape - a cohesive trainers tape developed for use in combination with PowerFlex. Together, PowerFlex and PowerTape create the Ultimate Taping System... the strongest taping system in existence.

The PowerFlex Taping System
JUST GOT BETTER.



TABLE OF CONTENTS:

Our Philosophy and History.....1

The Ultimate Taping System.....2

CREATIVE SOLUTIONS IN TAPING METHODS TO PREVENT INJURY

Ankle.....3-7

Upper Extremity.....8-10

Arch, Mid-Foot and Great Toe... 11-13

Ordering Information.....14

THE ULTIMATE TAPING SYSTEM



PowerTape™

Use in combination with PowerFlex.



THE ULTIMATE TAPING SYSTEM:

STEP 1:

Apply PowerFlex directly to the skin with proper un-wind tension.



STEP 2:

Cover PowerFlex with PowerTape.



- ✓ **PowerTape** is to be used in combination with PowerFlex.
- ✓ **PowerTape** replaces *traditional* cloth trainers tape.
- ✓ **PowerTape** combined with PowerFlex provides athletes with the STRONGEST taping system in existence.
- ✓ **PowerTape** contains less than .5% stretch - eliminating the "give" that exists when using *traditional* cloth trainers tape.
- ✓ **PowerTape** remains stable and comfortable providing maximum protection for the entire duration of games and practices.
- ✓ **PowerTape** tears by hand and sticks to itself for quick wrapping and superior staying power.
- ✓ **PowerTape** is sweat resistant and breathable. Taped areas will not slip.



2 - PATENTS PENDING

PowerTape™
NEW. INNOVATIVE. DIFFERENT.

TAPING METHODS TO PREVENT INJURY



PowerTape™

NEW. INNOVATIVE. DIFFERENT.



TAPING METHODS TO PREVENT INJURY

ANKLE - Preventive and Post-Injury Taping

The taping methods demonstrated here can be used for basic ankle injury prevention, ankle post-injury stability, and various post-injury situations. These techniques show multiple ways to provide injury prevention and post-injury stability for inversions, evasions and high ankle sprains.

I. REGULAR OR BASIC PREVENTIVE ANKLE



Figure 1.1



Figure 1.2



Figure 1.3



Figure 1.4

“PowerTape is the next natural step in providing your athletes with the most complete protective taping system available.”

Jerry Weber, A.T.C., P.T.
Head Athletic Trainer
University of
Nebraska - Lincoln

2.75" or 2" PowerFlex and 1.5" or 2" PowerTape. (See Figures 1.1 - 1.4)



Proudly Made in U.S.A.



PowerTape™

Use in combination with PowerFlex.



ANKLE - Preventive and Post-Injury Taping - Continued

II. PREVENTIVE ANKLE VARIATION - SPEED TAPING (Full)



Figure 1.5



Figure 1.6



Figure 1.7



Figure 1.8

2.75" or 2" PowerFlex and 1.5" or 2" PowerTape. (See Figures 1.5 - 1.8)



2 - PATENTS PENDING

PowerTape™
NEW. INNOVATIVE. DIFFERENT.

TAPING METHODS TO PREVENT INJURY



PowerTape™

NEW. INNOVATIVE. DIFFERENT.



TAPING METHODS TO PREVENT INJURY

ANKLE - Preventive and Post-Injury Taping - Continued

III. PREVENTIVE ANKLE VARIATION - SPEED TAPING (Modified)



Figure 1.9



Figure 1.10



Figure 1.11



Figure 1.12

“PowerFlex and PowerTape together allow the athletic trainer to apply a preventive and protective tape job that maintains it’s integrity longer and doesn’t get “soggy” with sweat, and thus is less likely to loosen.”

Jerry Weber, A.T.C., P.T.
Head Athletic Trainer,
University of Nebraska - Lincoln

2.75" or 2" PowerFlex and 1.5" or 2" PowerTape. (See Figures 1.9 - 1.12)





PowerTape™

Use in combination with PowerFlex.



ANKLE - Preventive and Post-Injury Taping - Continued

IV. PREVENTIVE ANKLE VARIATION - SPATTING



Figure 1.13



Figure 1.14



Figure 1.15



Figure 1.16

2.75" or 3" PowerFlex and 1.5" or 2" PowerTape. (See Figures 1.13 - 1.16)

PowerTape™

NEW. INNOVATIVE. DIFFERENT.



2 - PATENTS PENDING

TAPING METHODS TO PREVENT INJURY



PowerTape™

NEW. INNOVATIVE. DIFFERENT.



TAPING METHODS TO PREVENT INJURY

ANKLE - Preventive and Post-Injury Taping - Continued

V. PREVENTIVE ANKLE VARIATION - SPATTING (Modified/Speed Spatting)



Figure 1.17



Figure 1.18



Figure 1.19



Figure 1.20

2.75" or 3" PowerFlex and 1.5" or 2" PowerTape. (See Figures 1.17 - 1.20)

"I would definitely recommend trying PowerTape."

Andy Clawson, M.S., A.T.C., Director of Sports Medicine, *The Citadel*



Proudly Made in U.S.A.



PowerTape™

Use in combination with PowerFlex.

UPPER EXTREMITY - Preventive and Post-Injury Taping

The taping methods shown here can be used for wrist, hand, thumb, and elbow injury prevention, post injury stability and various post injury situations requiring special taping. These techniques provide post injury support for sprains and strains of these upper extremity areas.

I. REGULAR OR BASIC PREVENTIVE WRIST/HAND/THUMB



Figure 2.1



Figure 2.2



Figure 2.3

2" or 1.5" PowerFlex and 1.5" & 1" PowerTape. (Figure 2.1 - 2.3)

II. POST-INJURY WRIST VARIATION



Figure 2.4



Figure 2.5



Figure 2.6

2" or 1.5" PowerFlex and 1.5" & 1" PowerTape. (Figure 2.4 - 2.6)

PowerTape™
NEW. INNOVATIVE. DIFFERENT.



2 - PATENTS PENDING



PowerTape™

NEW. INNOVATIVE. DIFFERENT.



TAPING METHODS TO PREVENT INJURY

UPPER EXTREMITY - Preventive and Post-Injury Taping - Continued

III. REGULAR OR BASIC PREVENTIVE THUMB



Figure 2.7



Figure 2.8



Figure 2.9

2" or 1.5" PowerFlex, 1.5" & 1" PowerTape. (See Figures 2.7 - 2.9)

IV. POST-INJURY THUMB



Figure 2.10



Figure 2.11



Figure 2.12

2" or 1.5" PowerFlex, 1.5" & 1" PowerTape. (See Figures 2.10 - 2.12)

*"What sets PowerTape apart is that it contains virtually no stretch - less than .5%.
Traditional cloth trainers tape contains 5 - 6% stretch."*

Ron O'Neil, B.S., A.T.C.
Education, Research and Development
Andover Coated Products, Inc.



Proudly Made in U.S.A.



PowerTape™

Use in combination with PowerFlex.



UPPER EXTREMITY - Preventive and Post-Injury Taping - Continued

V. POST-INJURY ELBOW



Figure 2.13



Figure 2.14



Figure 2.15

2.75" or 3" PowerFlex and 1.5" or 2" PowerTape.
(Figure 2.13 - 2.15)

"Strongest post-injury taping system that provides maximum stability and range of motion control which is vital to preventing re-injury"

Ron O'Neil, B.S., A.T.C.
Education, Research and Development
Andover Coated Products, Inc.



2 - PATENTS PENDING

PowerTape™
NEW. INNOVATIVE. DIFFERENT.

TAPING METHODS TO PREVENT INJURY



PowerTape™

NEW. INNOVATIVE. DIFFERENT.



TAPING METHODS TO PREVENT INJURY

ARCH, MID-FOOT AND GREAT TOE - Post-Injury Taping

The taping methods here can be used for basic arch, mid-foot, and great toe injury prevention and post injury stability requiring special taping. These techniques provide post injury support for sprains and strains of the foot.

I. PREVENTIVE VARIATION - FOOT ARCH POST-INJURY



Figure 3.1



Figure 3.2



Figure 3.3

2.75" or 2" PowerFlex and
1.5" or 2" PowerTape.
(See Figures 3.1 - 3.3)

"The combination of PowerFlex and PowerTape has some unique characteristics that are very advantageous in the protective strapping of our athletes."

Andy Clawson, M.S., A.T.C.
Director of Sports Medicine
The Citadel



Proudly Made in U.S.A.



PowerTape™

Use in combination with PowerFlex.



ARCH, MID-FOOT AND GREAT TOE - Post-Injury Taping - Continued

II. POST-INJURY TURF TOE



Figure 3.4



Figure 3.5



Figure 3.6



Figure 3.7

2" PowerFlex, 1" PowerFlex and 1" or 1.5" PowerTape. (See Figures 3.4 - 3.7)

PowerTape™

NEW. INNOVATIVE. DIFFERENT.



2 - PATENTS PENDING

TAPING METHODS TO PREVENT INJURY



PowerTape™

NEW. INNOVATIVE. DIFFERENT.



ARCH, MID-FOOT AND GREAT TOE - Post-Injury Taping - Continued

III. POST-INJURY VARIATION - TURF TOE/ANKLE COMBINATION



Figure 3.8



Figure 3.9



Figure 3.10



Figure 3.11

“With so many positive performance and protection factors, it just makes sense to use the Andover combination of PowerFlex and PowerTape.”

Jerry Weber, A.T.C., P.T.
Head Athletic Trainer
University of
Nebraska - Lincoln

2.75" or 2" PowerFlex, 1" PowerFlex and 1", 1.5" or 2" PowerTape.
(See Figures 3.8 - 3.11)

TAPING METHODS TO PREVENT INJURY



ORDERING INFORMATION



Color Available:
White (WH)

PowerTape™ 2 - PATENTS PENDING

To be used in combination with PowerFlex for the Ultimate Taping System.

PowerTape replaces traditional cloth trainers tape and when used in combination with PowerFlex creates the strongest taping system in existence. PowerTape contains less than .5% stretch - eliminating the "give" that exists when using traditional trainers tapes.

CATALOG #	SIZE	ROLLS/CASE
ACP130-1"	1" x 15 yds.	48
ACP130-1.5"	1.5" x 15 yds.	32
ACP130-2"	2" x 15 yds.	24

PowerFlex U.S. PATENT NO. 5.762.623

PowerFlex is a patented cohesive flexible bandage that tears easily by hand for fast and efficient tape jobs. PowerFlex will not slip from sweat - it will remain in place. PowerFlex is lightweight and comfortable, yet provides the athlete with 23 lbs. of tensile strength for ultimate support. Available in 16 colors in sizes 1" to 6".



Colors Available:

Red (RD), Blue (BL), Green (GR),
Yellow (YL), White (WH), Black (BK),
Purple (PU), Neon Green (NG),
Neon Pink (NP), Light Blue (LB), Navy (NV),
Tan (TN), Teal (TE), Maroon (MR),
Orange (OR) and Gray (GY).

CATALOG #	SIZE	ROLLS/CASE
3710XX-048	1" X 6 yds.	48
3715XX-032	1.5" x 6 yds.	32
3720XX-024	2" x 6 yds.	24
3725XX-016	2.75" x 6 yds.	16
4725XX-048*	2.75" x 6 yds.	48
4730XX-048*	3" x 6 yds.	48
3730XX-016	3" x 6 yds.	16
3740XX-012	4" x 6 yds.	12
3760XX-008	6" x 6 yds.	8

When placing an order, please state catalog# and color.

Example: 3710WH. XX= Color Code.

* Player's Pack

Special Thanks

Andover would like to extend a special thanks to the following athletic training programs for their commitment to evaluating new products with an open mind. Their progressive thinking and methods have led to the development and refinement of PowerTape.

Jerry Weber and the athletic training staff at University of Nebraska.

Andy Clawson and the athletic training staff at The Citadel.



Andover Coated Products, Inc.
9 Fanaras Drive
Salisbury, MA 01952
Phone: 978.465.0044
Fax: 978.462.0003
Toll Free: 800.432.6686

For a distributor near you, call Toll Free: **800.432.6686**
Or, visit us on the web at: www.andovercoated.com

THE ULTIMATE TAPING SYSTEM

PowerTape™

NEW. INNOVATIVE. DIFFERENT.



ANDOVER COATED PRODUCTS, INC.
9 FANARAS DRIVE, SALISBURY, MA 01952
PHONE: 978.465.0044
TOLL FREE: 800.432.6686
FAX: 978.462.0003
www.andovercoated.com